



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Pararosaniline Acetate
Catalog# : S10143
CAS : 6035-94-5
CI number : 42500
Formula : $C_{21}H_{21}N_3O_2$
MW : 347.42

TEST	SPECIFICATION
Appearance	Dark Maroon Powder
Melting Point	203 ° C
Lambda Max	545 nm in 50 % Ethanol
Dye Content	85 % minimum

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

