



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Pararosaniline HCL Certified  
Catalog# : S10144  
CAS : 569-61-9  
CI number : 42500  
Formula :  $C_{19}H_{17}N_3 \cdot HCl$   
MW : 323.82

TEST	SPECIFICATION
Appearance	Green crystalline powder
Assay	88 % maximum
Melting Point	268 - 270 ° C
Solubility: 1 mg/ml, EtOH	Clear, dark red solution
Lambda max	544nm in 50% ethanol

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

