



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Safranin O Certified
Catalog# : S10174
CAS : 477-73-6
CI number : 50240
Formula : $C_{20}H_{19}ClN_4$
MW : 350.85

TEST	SPECIFICATION
Appearance	Dark-Green crystalline Powder
Lambda max	530 nm in 50 % Ethanol
Dye content	95% minimum
Solubility: 1 mg/ml, H ₂ O	Opaque red solution
Loss on drying	5% maximum

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

