



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Sudan II
Catalog# : S10178
CAS : 3118-97-6
CI number : 12140
Formula : $C_{18}H_{16}N_2O$
MW : 276.34

TEST	SPECIFICATION
Appearance	Orange Powder
Lambda max	493 nm
Melting Point	156 - 158 °C
Dye Content	85 % minimum

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

