



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Sudan III Certified  
Catalog# : S10179  
CAS : 85-86-9  
CI number : 26100  
Formula : C<sub>22</sub>H<sub>16</sub>N<sub>4</sub>O  
MW : 352.40

TEST	SPECIFICATION
Assay	95% minimum
Appearance	Red-Brown Powder
Lambda max	507 (354) nm (in Toluene)
Melting Point	199 ° C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

