



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Sudan IV Certified
Catalog# : S10180
CAS : 85-83-6
CI number : 26105
Formula : C₂₄H₂₀N₄O
MW : 380.45

TEST	SPECIFICATION
Description	Dark Red-Brown Powder
Identification	Conforms to structure and standard
Dye Content	> 80.0%
Melting Point	199 C (dec.)
Spectroscopy	Peak 520 nm in Toluene
Solubility (1 mg/ml, CHCL ₃)	Opaque dark red solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

