



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Sudan Black B Certified
Catalog# : S10181
CAS : 4197-25-5
CI number : 26150
Formula : $C_{29}H_{24}N_6$
MW : 456.55

TEST	SPECIFICATION
Appearance	Dark brownish black powder
UV absorption @ 596 nm	≥ 20000
Loss on drying	$\leq 5.0\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

