



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Thionin Acetate Certified
Catalog# : S10189
CAS : 78338-22-4
CI number : 52000
Formula : $C_{12}H_{10}N_3S.C_2H_3O_2$
MW : 287.33

TEST	SPECIFICATION
Appearance	Very Dark Green to Bluish Black Powder
Lambda max	598 nm (In Water)
Dye Content	85% minimum
Loss on drying	7% maximum
To be certified by Biological Stains Commission	

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

