



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Wrights Stain Certified  
Catalog# : S10197  
CAS : 68988-92-1  
CI number : -  
Formula : -  
MW : -

TEST	SPECIFICATION
Description	Green Crystalline Powder
Solubility 0.1% (MeOH)	Conforms Opaque c= 1mg/ml: Methanol
Solubility (Color)	Dark Blue
Infrared spectrum	Conforms to structure
Loss on drying At 105 deg. C/1hour	
HPLC:	
Suitable as a Biological Stain:	Passes Test
Certified by Biological Stain Commission	Passes Test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

