



Product Specification

Product Information:

Products:	Crystal violet staining solution
Catalog:	SS10270
Grade:	Stain solution

Reagent Used:

Reagents	CAS #	Grade
Ethyl alcohol	64-17-5	ACS
Crystal violet	548-62-9	ACS
Oxalic acid ammonium salt monohydrate	6009-70-7	ACS

Specification:

Test	Specification
Appearance	Purple to very dark purple solution
Suitability for use	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

