



Product Specification

Product Information:

Products:	Eosin Y (0.5% Aqueous solution)
Catalog:	SS10272
Grade:	Stain solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Ethyl alcohol	64-17-5	ACS
Eosin Y	17372-87-1	ACS

Specification:

Test	Specification
Appearance	Dark red color liquid
Suitability for Use	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

