



Product Specification

Product Information:

Products:	Gram's Iodine staining
Catalog:	SS10277
Grade:	Stain solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Potassium iodide	7681-11-0	ACS
Iodine	7553-56-2	ACS

Specification:

Test	Specification
Appearance	Dark brown liquid
Suitability for use	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

