



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	Methylene blue (1% Aqueous solution)
Catalog:	SS10293
Grade:	Stain solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Methylene blue trihydrate	7220-79-3	ACS

Specification:

Test	Specification
Appearance	Blue liquid.
UV/Visible spectrum	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

