



Tri Stains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

### Product Information:

<b>Products:</b>	Safranine solution
<b>Catalog:</b>	SS10307
<b>Grade:</b>	Stain solution

### Reagent Used:

Reagents	CAS #	Grade
Ethyl alcohol	64-17-5	ACS
Methyl alcohol	67-56-1	ACS
Safranine	477-73-6	Certified
Water	7732-18-5	USP/ASTM

### Specification:

Test	Specification
Appearance	Dark red liquid
Suitability for use	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

