



## Product Specification

### Product Information:

<b>Products:</b>	Lactophenol cotton blue
<b>Catalog:</b>	SS10322
<b>Grade:</b>	Stain solution

### Reagent Used:

Reagents	CAS #	Grade
Phenol	108-95-2	ACS
Glycerin	56-81-5	ACS
Lactic acid	79-33-4	ACS
Water	7732-18-5	USP/ASTM
Cotton blue	28983-56-4	Certified

### Specification:

Test	Specification
Appearance	Blue to very dark blue liquid
Suitability for Use	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

