



## Product Specification

### Product Information:

<b>Products:</b>	Potassium hydroxide 3%
<b>Catalog:</b>	SS10492
<b>Grade:</b>	Stain solution

### Reagent Used:

Reagents	CAS #	Grade
Potassium hydroxide	1310-58-3	ACS
Water	7732-18-5	USP/ASTM

### Specification:

Test	Specification
Appearance	Colorless liquid
Assay	3.0%±1.0%
Suitability for Use	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

