

# Product Specification

**Product Information:**

<b>Products:</b>	Acetic acid Ammonium acetate buffer TS
<b>Catalog:</b>	TS14004
<b>Grade:</b>	Test solution

**Reagent Used:**

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Ammonium acetate	631-61-8	ACS
Acetic acid glacial	64-19-7	ACS

**Specification:**

Test	Specification
Appearance	Colorless liquid.
pH at 25°C	4.5 to 4.9

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

