

Product Specification

Product Information:

Products:	Alkaline Cupric Tartrate TS
Catalog:	TS14019
Grade:	Test solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Copper sulfate pentahydrate	7758-99-8	ACS
Potassium sodium tartrate tetrahydrate	6381-59-5	ACS
Sodium hydroxide	1310-73-2	ACS

Specification:

Test	Specification
Appearance	Clear blue color liquid
Suitability for use	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

