



Product Specification

Product Information:

Products:	Ammonia-Cyanide TS
Catalog:	TS14030
Grade:	Test solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Ammonium hydroxide	1336-21-6	ACS
Potassium Cyanide	151-50-8	ACS

Specification:

Test	Specification
Appearance	Colorless liquid.

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

