

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	Ammonium carbonate TS
Catalog:	TS14034
Grade:	Test solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Ammonium hydroxide	1336-21-6	ACS
Ammonium carbonate	506-87-6	ACS

Specification:

Test	Specification
Appearance	Colorless liquid.
Assay	4.73 to 5.19 N at 20°C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

