

# Product Specification

## Product Information:

<b>Products:</b>	Ammonium Molybdate TS
<b>Catalog:</b>	TS14038
<b>Grade:</b>	Reagent solution

## Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Nitric Acid	7697-37-2	ACS
Molybdic Acid	7782-91-4	ACS
Ammonium Hydroxide	1336-21-6	ACS

## Specification:

Test	Specification
Appearance	Colorless to pale yellow-green liquid
Suitability for use	Yellow precipitate with Phosphate

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

