

Product Specification

Product Information:

Products:	Ammonium oxalate TS
Catalog:	TS14039
Grade:	Reagent solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Ammonium oxalate	6009-70-7	ACS

Specification:

Test	Specification
Appearance	Colorless liquid
Assay	3.47-3.53 % (w/v)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

