

# Product Specification

## Product Information:

<b>Products:</b>	Barium chloride TS
<b>Catalog:</b>	TS14049
<b>Grade:</b>	Test solution

## Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Barium chloride dihydrate	10326-27-9	ACS

## Specification:

Test	Specification
Appearance	Colorless liquid.
Assay	11.7 to 12.3 % (w/v)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

