

Product Specification

Product Information:

Products:	Bromocresol green TS
Catalog:	TS14061
Grade:	Reagent Solution

Reagent Used:

Reagents	CAS #	Grade
Ethyl Alcohol	64-17-5	ACS
Methyl Alcohol	67-56-1	ACS
Bromocresol green	76-60-8	ACS

Specification:

Test	Specification
Appearance	Dark blue liquid
Suitability for Use	Blue (pH > 5.4) – Green (pH 4.5) – Yellow (pH < 3.8)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

