

Product Specification

Product Information:

Products:	m-Cresol Purple TS
Catalog:	TS14077
Grade:	Test solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
m-Cresol purple sodium salt	62625-31-4	ACS

Specification:

Test	Specification
Appearance	Reddish orange
pH at 25°C	6.85 – 7.05
Suitability for use	Purple (pH > 9.0) – Yellow (pH 2.8 to pH 7.4) – Red (pH < 1.2)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

