

Product Specification

Product Information:

Products:	Cupric Citrate TS
Catalog:	TS14084
Grade:	Reagent Solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Sodium carbonate	497-19-8	ACS
Citric acid	77-92-9	ACS
Copper sulfate pentahydrate	7758-99-8	ACS

Specification:

Test	Specification
Appearance	Blue liquid
Refractive index	To Pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

