

Product Specification

Product Information:

Products:	Erichrome Black TS
Catalog:	TS14107
Grade:	Test solution

Reagent Used:

Reagents	CAS #	Grade
Methyl alcohol	67-65-1	ACS
Hydroxylamine hydrochloride	5470-11-1	ACS
Eriochrome black T	1787-61-7	ACS

Specification:

Test	Specification
Appearance	Deep red liquid.
Suitability for use	Wine red (Excess Ca/Mg) – Blue (Excess EDTA)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

