



# Product Specification

**Product Information:**

<b>Products:</b>	Ferric ammonium sulfate TS
<b>Catalog:</b>	TS14110
<b>Grade:</b>	Test solution

**Reagent Used:**

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Ferric ammonium sulfate dodecahydrate	7783-83-7	ACS

**Specification:**

Test	Specification
Appearance	Amber liquid
Suitability for use	7.9 to 8.1 % (w/v)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

