

Product Specification

Product Information:

Products:	Iodobromide TS
Catalog:	TS14137
Grade:	Reagent Solution

Reagent Used:

Reagents	CAS #	Grade
Acetic acid	64-19-7	ACS
Iodine monobromide	7789-33-5	Practical

Specification:

Test	Specification
Appearance	Deep red-orange liquid
Total Halogen content	0.191 – 0.195 N at 20°C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

