

Product Specification

Product Information:

Products:	Lead Acetate TS
Catalog:	TS14143
Grade:	Reagent Solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Lead acetate trihydrate	6080-56-4	ACS

Specification:

Test	Specification
Appearance	Colorless liquid
Assay	9.46 – 9.54 % (w/v)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

