



Product Specification

Product Information:

Products:	Methyl orange TS
Catalog:	TS14167
Grade:	Reagent Solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Methyl Orange	547-58-0	ACS

Specification:

Test	Specification
Appearance	Orange liquid
Suitability for Use	Yellow (pH > 4.4) – Pinkish-Orange (pH 3.2) – Pink (pH < 3.0)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

