

Product Specification

Product Information:

Products:	Methyl red TS 2
Catalog:	TS14170
Grade:	Test solution

Reagent Used:

Reagents	CAS #	Grade
Methyl alcohol	67-56-1	ACS
Ethyl alcohol	64-17-5	ACS/Reagent
Isopropyl alcohol	67-63-0	ACS
Water	7732-18-5	ASTM/USP
Methyl red	493-52-7	ACS
Sodium hydroxide	1310-73-2	ACS

Specification:

Test	Specification
Appearance	Dark red liquid
Suitability for use	Pink (pH <4.2) -Yellow (pH > 6.2)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

