

Product Specification

Product Information:

| | |
|------------------|------------------|
| Products: | Ninhydrin TS |
| Catalog: | TS14189 |
| Grade: | Reagent solution |

Reagent Used:

| Reagents | CAS # | Grade |
|-----------|-----------|----------|
| Water | 7732-18-5 | ASTM/USP |
| Ninhydrin | 485-47-2 | ACS |

Specification:

| Test | Specification |
|---------------------|--|
| Appearance | Colorless to light yellow liquid. |
| Suitability for use | Blue color in presence of amino acids, proteins and peptides |

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

