

Product Specification

Product Information:

Products:	Orthophenanthroline TS
Catalog:	TS14193
Grade:	Test Solutions

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Ferrous sulfate heptahydrate	7782-63-0	ACS
1,10-Phenanthroline monohydrate	5144-89-8	ACS

Specification:

Test	Specification
Appearance	Dark red liquid
Suitability for Use	Green to Bluish Green (Reduced)- Redish (Oxidized)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

