

Product Specification

Product Information:

| | |
|------------------|----------------------|
| Products: | pH 4.7 Phenol red TS |
| Catalog: | TS14199 |
| Grade: | Test solution |

Reagent Used:

| Reagents | CAS # | Grade |
|------------------|-----------|----------|
| Water | 7732-18-5 | ASTM/USP |
| Acetic acid | 64-19-7 | ACS |
| Sodium hydroxide | 1310-73-2 | ACS |
| Phenol red | 143-74-8 | ACS |
| Ammonium sulfate | 7783-20-2 | ACS |

Specification:

| Test | Specification |
|--------------------------|------------------|
| Appearance | Yellow liquid. |
| pH at 25°C | 4.65 to 4.75 |
| Sensitivity as Indicator | To Pass the test |

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

