

Product Specification

Product Information:

Products:	Potassium Acetate TS
Catalog:	TS14212
Grade:	Reagent solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Potassium Acetate	127-08-2	ACS

Specification:

Test	Specification
Appearance	Colorless liquid
Refractive Index	To Pass Test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

