



# Product Specification

**Product Information:**

<b>Products:</b>	Potassium Iodide TS
<b>Catalog:</b>	TS14222
<b>Grade:</b>	Test solution

**Reagent Used:**

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Potassium Iodide	7681-11-0	ACS

**Specification:**

Test	Specification
Appearance	Colorless to light yellow liquid
Assay	16.4 to 16.6% (w/v)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

