

Product Specification

Product Information:

Products:	Potassium Pyroantimonate TS
Catalog:	TS14226
Grade:	Test solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Potassium hydroxide	1310-58-3	ACS
Potassium pyroantimonate	12208-13-8	Reagent
Sodium hydroxide	1310-73-2	Reagent

Specification:

Test	Specification
Appearance	Colorless liquid.
Suitability for use	Dense precipitate with sodium

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

