

Product Specification

Product Information:

Products:	Silver diethyldithiocarbamate TS
Catalog:	TS14240
Grade:	Test solution

Reagent Used:

Reagents	CAS #	Grade
Silver diethyldithiocarbamate	1470-61-7	ACS
Pyridine	110-86-1	ACS

Specification:

Test	Specification
Appearance	Intensely yellow color liquid.
Suitability for use	To pass the test
Assay	

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

