



Product Specification

Product Information:

Products:	Sodium hypochlorite TS
Catalog:	TS14260
Grade:	Reagent Solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Sodium hypochlorite	7681-52-9	Commercial

Specification:

Test	Specification
Appearance	Colorless to greenish-yellow liquid
Assay	4.0 to 6.0%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

