



Product Specification

Product Information:

Products:	Starch TS (Iodide Free)
Catalog:	TS14274
Grade:	Reagent Solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Starch soluble	9005-84-9	ACS

Specification:

Test	Specification
Appearance	White translucent liquid.
Suitability for use	Colorless (Iodine absent) – Blue (Iodine present)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

