

# Product Specification

**Product Information:**

<b>Products:</b>	Fehling solution A TS
<b>Catalog:</b>	TS14412
<b>Grade:</b>	Test solution

**Reagent Used:**

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Copper sulfate pentahydrate	7758-99-8	ACS

**Specification:**

Test	Specification
Appearance	Blue liquid
Assay	0.2770 to 0.2780 N at 20°C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

