



# Product Specification

**Product Information:**

<b>Products:</b>	Fehling solution B TS
<b>Catalog:</b>	TS14413
<b>Grade:</b>	Test solution

**Reagent Used:**

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Potassium sodium tartrate tetrahydrate	6381-59-5	ACS
Sodium hydroxide	1310-73-2	

**Specification:**

Test	Specification
Appearance	Colorless liquid

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

