

Product Specification

Product Information:

Products:	Ferroxyl test solution Kit
Catalog:	TS14501
Grade:	Reagent solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Nitric Acid	7697-37-2	ACS
Potassium ferricyanide	13746-66-2	ACS

Specification:

Test	Specification
Appearance	Slight yellow to pale yellow-green liquid
Suitability for use	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

