



## Product Specification

**Product Information:**

<b>Products:</b>	Cupriethylenediamine 1.0M
<b>Catalog:</b>	VS17042
<b>Grade:</b>	Standard volumetric solution

**Reagent Used:**

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Cupriethylenediamine	13426-91-0	ACS

**Specification:**

Test	Specification
Appearance	Dark blue
Assay - Molarity	0.98 – 1.02M
Ratio Cu: Ethylenediamine	1.95 – 2.05
Viscosity	1.147 – 1.167 cSt

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

