



Product Specification

Product Information:

Products:	EDTA solution 0.05M
Catalog:	VS17045
Grade:	Standard volumetric solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Ethylenediaminetetraacetic acid (EDTA) disodium salt	6381-92-6	ACS

Specification:

Test	Specification
Clarity @25°C	Clear solution
Assay	0.05 M \pm 5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

