

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	Hydrochloric acid 0.5N
Catalog:	VS17064
Grade:	Standard volumetric solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Hydrochloric acid	7647-01-0	ACS

Specification:

Test	Specification
Clarity @25°C	Clear solution
Assay	0.5 N \pm 5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

