



## Product Specification

**Product Information:**

<b>Products:</b>	Perchloric acid 0.1N in Acetic acid
<b>Catalog:</b>	VS17106
<b>Grade:</b>	Standard volumetric solution

**Reagent Used:**

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Acetic Acid	64-19-7	ACS
Acetic Anhydride	108-24-7	ACS
Perchloric Acid	7601-90-3	ACS

**Specification:**

Test	Specification
Appearance	Clear Colorless to yellow solution
Water	0.02% to 0.5%
Assay	0.1 N $\pm$ 5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

