



## Product Specification

**Product Information:****Products:** Potassium hydroxide 0.1N (0.1M)**Catalog:** VS17138**Grade:** Standard volumetric solution**Reagent Used:**

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Potassium hydroxide	1310-58-3	ACS

**Specification:**

Test	Specification
Clarity @25°C	Clear solution
Assay	0.1 N $\pm$ 5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

