

Product Specification

Product Information:

Products:	Sodium thiosulfate 0.05N
Catalog:	VS17250
Grade:	Standard volumetric solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Sodium thiosulfate pentahydrate	10102-17-7	ACS
Sodium carbonate	497-19-8	ACS

Specification:

Test	Specification
Appearance	Clear solution
Assay at 20°C	0.05 N ± 5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

